

## Dia-Temp™

Nano-silver contained light-curing temporary filling material

**Dia-Temp™** is a light curing temporary filling material and it includes nano-sized silver.

**Dia-Temp™** is for the temporary restoration of inlay and onlay preparations as well as the prefabricated temporary crowns and bridges. With high elasticity and surface hardness, **Dia-Temp™** demonstrates also good adhesion.

**Dia-Temp™** is used without any additional mixing and also removed all at once without any particles in the cavity.

### Composition

- Polyurethane dimethacrylate
- Hydrophilic methacrylate
- Nano-sized silica
- Nano-sized silver
- Catalysts and stabilizer

### Indications

- Temporary restoration of inlay and onlay preparations
- Temporary fillings of all kinds
- Temporary fillings and sealings

### Recommended Procedure

1. Clean the teeth with water gently and isolate the teeth with a rubber dam.
2. Remove all decay from the tooth and prepare the cavity.
3. Dispense necessary amount of **Dia-Temp™** with a dental instrument such as a spatula and etc., and fill it in the cavity.
4. Light cure **Dia-Temp™** for 40 seconds with a standard light curing unit. (DiaDent **D-Lux** is recommended.) If light curing unit output is below 400mW/cm<sup>2</sup> by radiometer, increased curing time is needed.
5. Check the occlusion and remove excess of **Dia-Temp™** on the surface.
6. Remove excess of **Dia-Temp™** on the surface with an explorer or reasonable dental instrument.

### Package

**Dia-Temp™** kit (3g syringe × 3)  
Blue (2003-1101) Yellow (2003-1201)

### Contact Information

#### DiaDent Group International

626, Yeonje-ri, Gangoe-myeon, Cheongwon-gun,  
Chungcheongbuk-do, Korea. 363-951  
Tel : 82-43-266-2315  
Fax : 82-43-262-8658  
http://www.diadent.co.kr  
E-mail : diadent@chollian.net

#### DiaDent Group International Inc.

11-3871 North Fraser Way, Burnaby  
B.C., Canada, V5J 5G6  
Tel : 1-604-451-8851  
Fax : 1-604-451-8865  
http://www.diadent.com  
E-mail : diadent@diadent.com  
Toll Free : 1-877-DIA-DENT(1-877-342-3368)

#### DiaDent Europe B.V.

Antennestraat 70, 1322AS  
Almere, the Netherlands  
Tel : 31-36-549-8607  
Fax : 31-36-536-7317  
http://www.diadenteurope.com  
E-mail : diadent@diadenteurope.com

Copyright Mar.21, 2012 [IFU: Diapex plus(rev.1)]

### Precautions

- **Dia-Temp™** is used by a dentist and dental professional for dental purpose only and Instruction for Use is carefully read prior to using the **Dia-Temp™**.
- To use cold **Dia-Temp™** from the refrigerator, please keep it approx. 30 minutes at room temperature and then apply to use the **Dia-Temp™**. Close the syringe cap quickly and tightly after cutting out necessary amount of **Dia-Temp™**.
- Do not use after expiry date.
- In case of slow or non-polymerization, try to check the light curing unit firstly, and then contact the nearest dealer(s) or manufacturer directly regarding same quality matters.
- Using protective gloves is recommended.
- In accidental contact with eyes or prolonged contact with oral tissue occurs, flush immediately with a large amount of water.

### Symbols




Use by

LOT

Batch Code



Consult Instructions for use



Temperature  
Limitation



Manufacturer

EC REP

Authorised Representative  
in the European Community

### Warranty

DiaDent Group International will replace product that is proven to be defective. DiaDent Group International does not accept liability for any damage or loss, direct or consequential, stemming from the use of or inability to use the product described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.