

 **semper**
med

derma plus



Produced for increased safety

Manufactured while using the specially developed double-dipping process, the sempermed® derma plus provides excellent donning properties.



High tear resistance



Secure grip



**Anatomically
perfect shape**



derma plus



High tear resistance

Better protection –

The typical force at break of the glove is significantly higher than the standard requirement and therefore provides the user better protection.



Secure grip

Through micro-rough surface –

Good grip through micro-rough surface texture, even in a moist environment.





Anatomically perfect shape

Comfortable fit –

With broader back of the hand and naturally curved fingers.

Technical data

Type	sterile surgical glove for single use, lightly powdered	
Color	natural white	
Size / Overall length as per EN 455-2	5½, 6 and 6½ / 7, 7½ and 8 / 8½ and 9	median 270 mm median 280 mm median 285 mm
Wall thickness measured in single layer in the palm area	0.19 - 0.24 mm	
Barrier performance	AQL 0.65	
Material	natural rubber latex (NR)	
Glove shape	fully anatomical with rolled rim	
Force at break as per EN 455-2	≥ 12 N, exceeding EN 455-2 (≥ 12 N instead of 9 N)	
Shelf life in original packaging if stored according to product specification	5 years	
Sterilization	radiation STERILE R	
Labeling & conformity to standards	EN 455-1/-2/-3/-4, EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, ISO 374-5, EN 421, EN 556, ISO 2859-1, ISO 15223-1, ASTM F1671/F1671M	
Purpose	<p>Single-use medical device class IIa according to MD Directive 93/42/EEC</p> <p>Single-use protective glove – PPE category III (protection against chemical substances for limited time)* according to PPE Regulation (EU) 2016/425</p>	 



Packaging

- Left and right glove with turned up cuff in inner pouch, ozone-tight, sealed in peel pack.
- In dispenser carton with sterilization indicator: 50 pairs
- In transport carton with sterilization indicator: 300 pairs

Article numbers

size 5½	823152529
size 6	823152609
size 6½	823152629
size 7	823152709
size 7½	823152729
size 8	823152809
size 8½	823152829
size 9	823152909