

Lilac Utility Gloves

For use in the sterilization area, and during operatory breakdown and disinfection processes.

- Nitrile, lilac utility gloves stand out to professionals as the color of dental hygiene
- Flocklined for increased ease of donning and doffing
- Resist chemicals better than exam gloves
- Puncture resistant for improved protection
- Molded texture improves tactile feel and grip
- Can be autoclaved up to 250° F/121° C up to five times

Not made with natural rubber latex



SIZE	PART CODES (three pairs per package)
Small (Size 7)	40-060
Medium (Size 8)	40-062
Large (Size 9)	40-064
ExtraLarge (Size 10)	40-066




The cleaning process is an important step in the Instrument Management Cycle. Make it safer and more efficient with the Instrument Management System (IMS).

The CDC recommends the use of puncture and chemical-resistant utility gloves when cleaning instruments and performing housekeeping tasks involving contact with blood, other potentially infectious materials, or chemicals.

USE THIS CHART TO DETERMINE YOUR APPROPRIATE GLOVE SIZE

Place hand as indicated to select correct size in order to avoid hand fatigue. It is important that each clinician have their own pair of well-fitting utility gloves.



For more information: Call 1-800-Hu-Friedy, contact your authorized HuFriedyGroup representative, or visit us online at Hu-Friedy.com

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC its affiliates or related companies, unless otherwise noted.

©2020 Hu-Friedy Mfg. Co., LLC. All rights reserved. IMS-340 / 1020