


M5™ Thermal
Copper Nickel Titanium

Precision engineered and manufactured by G&H® in the U.S.A.



ORTHODONTICS®
Quality Manufacturing. Exceptional Service.



GHOrthodontics.com

TOLL-FREE 800-526-1026


orders@ghortho.com

G&H Orthodontics 40 Linville Way Franklin, IN 46131 U.S.A. 800-526-1026 U.S. & Canada +1 317-346-6655 International


M5™ Thermal Copper Nickel Titanium

Precision engineered and manufactured by G&H® in the U.S.A.

M5™ Thermal Copper Nickel Titanium*



ANSI/ADA Specifications #32, Orthodontic Wires: 2017
Force values expressed at 3mm deflection.




G&H QUALITY MEANS:

- Gentler loading forces for easier engagement
- Precise temperature transformation of 35°C
- Consistent and predictable tooth moving forces
- Proprietary mechanical high polish technique for reduced surface friction
- US manufacturing quality



Call **800-526-1026** to order
or email us at orders@ghortho.com

Copper Nickel Titanium Wires
Force measured from 3-pt bend test



See and order our full line of products at GHOrthodontics.com